



Meshindia
Analytics Private Limited

PUREDROP - WATER PURIFICATION SYSTEMS

2025

The PUREDROP Series Offers A
Versatile Range Of Ultrapure
Water Systems Tailored For
Laboratory Applications.

www.meshindiaanalytics.com

Introduction



- ✓ The PUREDROP Series offers a versatile range of ultrapure water systems tailored for laboratory applications. Each model is designed to meet specific user needs, from small-volume HPLC use to multi-grade water production for regulated industries.
- ✓ The industry's most complete water purification solution for highly regulated applications.
- ✓ When your application requires the ultimate in water purity, PUREDROP provides the perfect solution. Consistently delivering water purity of 18.2 MΩ·cm and underpinned by advanced technologies, the PUREDROP enables you to focus on attaining accurate results while ensuring an uninterrupted workflow.

Typical Scientific Applications


- ICP-MS (Inductively Coupled Plasma Mass Spectrometry)
- Molecular biology techniques
- Ultra-trace analysis
- Electrochemistry
- Electrophoresis
- GFAAS (Graphite Furnace Atomic Absorption Spectrophotometry)
- HPLC
- IC (Ion Chromatography)
- ICP-AES (Inductively Coupled Plasma Atomic Emission Spectrometry)
- Mammalian and bacterial cell culture
- Molecular biology
- Plant tissue culture
- Qualitative analysis

PUREDROP
Mini



Features of Mini

- Many Laboratories require ultrapure water everyday- for HPLC or for other instrumental analysis- but only as small volumes at a time, up to a total of a few liters at the most.
- However, this is no reason to do without professional, freshly prepared ultrapure water and buy expensive HPLC water. After all with the repeated taking of small volumes from an HPLC water container, the quality of the water and the TOC-value no longer fulfill the specifications. The ultrapure water is no longer ultrapure!
- Always fresh. Always available in the volume required. Always with the required quality. And always for only a few cents per liter!
- Compact and efficient for labs requiring small volumes.
- Fresh ultrapure water on demand, ideal for HPLC.
- Internal 6L tank — no external tank needed.
- Cost-effective alternative to bottled HPLC water.



| | |
|------------------------------|-----------------------------------|
| 1. Solenoid Valve | 2. Conductivity Sensor |
| 3. Pre-treatment Module | 4. Boost Pump |
| 5. RO Module | 6. Ball Valve |
| 7. Dual Wave UV Cartridge | 8. Purification Module |
| 9. Ultra-Purification Module | 10. Micro-Filtration Cartridge |
| 11. Point-of-use Filter | 12. Internal Water Reservoir - 6L |

| Model | PUREDROP Mini |
|---|-------------------------------|
| Feed Water Requirement | |
| Source | Tap water |
| Conductivity* | <2000 µs/cm |
| Hardness** | <450 ppm as CaCO ₃ |
| Pressure | 0.05–0.5 MPa (7–72 psi) |
| Temperature | 5–40°C |
| Purification Water (Class III) | |
| Ionic Rejection | >95% |
| Bacteria Rejection | >99% |
| Conductivity | 1–20 µs/cm (RO2–5 µs/cm) |
| Productivity Rate | 30 L/h |
| Ultra-purification Water (Class I) | |
| Resistivity at 25°C | 18.2 MΩ.cm |
| Conductivity at 25°C | 0.055 µs/cm |
| TOC Level*** | 1–5 ppb |
| Endotoxin (Pyrogens) | N/A |
| Particulate (>0.02 µm) | <1 pc/ml |
| Bacteria*** | <0.1 cfu/ml |
| Rnase/Dnase** | N/A |
| Ultra-purification Water (Class I) | |
| Manual dispense flow rate | 1.2–1.5 L/min |
| Automatic dispense volume | 100–60000 ml |
| Electrical Requirements | |
| Electrical Voltage | 220V ±10% |
| Electrical Frequency | 50Hz |
| Net Weight - Main units | 32 kg |
| Net Weight - Water tank | N/A |
| External Dimensions (Main units) | 315 × 525 × 570 mm |
| External Dimensions (Water tank) | N/A |

- If feed water quality is poor (Conductivity >1000 µs/cm), strengthened pretreatment module and RO-2 type is highly recommended
- When hardness of feed water is high (>450 ppm as CaCO₃), 0.5T water softener is recommended
- Dual wave UV module need to be adopted. Also dependent on feed water, recommended feed TOC <30 ppb


PUREDROP-N



Features of PUREDROP-N

- The industry’s most complete water purification solution for highly regulated applications.
- When your application requires the ultimate in water purity, PUREDROP N Series provides the perfect solution. Consistently delivering water purity of 18.2MΩ.cm and underpinned by the advanced technologies, the PUREDROP N enables you to focus on attaining accurate results while ensuring an uninterrupted work flow.
- Designed for highly regulated environments.
- Dual-wave UV ensure <5ppb TOC.
- Remote dispensing to reduce contamination.

| | |
|-------------------------------|------------------------------|
| 1 Solenoid Valve | 2 Conductivity Sensor |
| 3 Pre-treatment Module | 4 Boost Pump |
| 5 RO Module | 6 Water Tank |
| 7 Ball Valve | 8 Dual Wave UV Cartridge |
| 9 Purification Module | 10 Ultra-Purification Module |
| 11 Ultra-Filtration Cartridge | 12 Point-Of-Use Filter |
| 13 Remote Water Dispenser | |



| Specification | Details |
|---|------------------------------|
| Model | PUREDROP N 15UV |
| Feed Water Requirement | |
| Source | Tap water |
| Conductivity* | <2000us/cm |
| Hardness** | <450ppm as CaCO ₃ |
| Pressure | 0.05~0.5MPa (7~72psi) |
| Temperature | 5~40°C |
| Purification Water (Class III) | |
| Ionic Rejection | >95% |
| Bacteria Rejection | >99% |
| Conductivity | 1~20us/cm (RO-2 5us/cm) |
| Productivity Rate | 15L/h |
| Ultra-purification Water (Class I) | |
| Resistivity At 25°C | 18.2MΩ.cm |
| Conductivity At 25°C | 0.055us/cm |
| TOC Level*** | 1~5ppb |
| Endotoxin (Pyrogens)**** | <0.001 EU/ml |
| Particulate (≥0.02um) | <1pc/ml |
| Bacteria*** | <0.1cfu/ml |
| Rnase / Dnase** | N/A |
| Manual dispense flow rate | 1.5~2.0L/min |
| Automatic dispense volume | 100~600000ml |
| Electrical Requirements | |
| Electrical Voltage | 110V/220V ±10% |
| Electrical Frequency | 50Hz/60Hz |
| Net Weight | |
| Main units | 24kg |
| Water tank (30L) | 7kg |
| External Dimensions (W×D×H) | |
| Main units | 315×525×570mm |
| Water tank (30L) | 380×380×595mm |

- If feed water quality is poor (Conductivity > 1000 μs/cm), strengthened pretreatment module and RO-2 type is highly recommended
- When hardness of feed water is high (> 450 ppm as CaCO₃), 0.5 T water softener is recommended
- Dual wave UV module needs to be adopted. Also dependent on feed water, recommended feed TOC < 30 ppb
- Ultra-filtration module needs to be adopted. Feed water needs to be satisfied as above


PUREDROP + - NE/NET



Features of PUREDROP + - NE/NET

- The industry’s most complete water purification solution for highly regulated applications.
- When your application requires the ultimate in water purity, PUREDROP-NE Series provides the perfect solution. Consistently delivering water purity of 18.2MΩ.cm and underpinned by the advanced technologies, the PUREDROP-NE enables you to focus on attaining accurate results while ensuring an uninterrupted work flow.
- Delivers Class I, II, and III water in one system.
- Class II high-purity stage for intermediate-grade use.
- Enhanced organic compound removal.
- Ideal for research centers with diverse purity needs.
- Dual-wave UV and TOC monitoring ensure <5ppb TOC.
- Includes EDI and ultrafiltration modules for maximum purity.

| | |
|-------------------------------|-------------------------------|
| 1 Solenoid Valve | 2 Conductivity Sensor |
| 3 Pre-treatment Module | 4 Boost Pump |
| 5 RO Module | 6 Water Tank |
| 7 Ball Valve | 8 Dual Wave UV Cartridge |
| 9 Ultra - Purification Module | 10 Ultra-Filtration Cartridge |
| 11 Point-Of-Use Filter | 12 Remote Water Dispenser |
| 13 TOC monitor cell* | 14 EDI Module |



PUREDROP + - NE/NET Specifications

Feed Water Requirement

| Parameter | Specification |
|---------------|-------------------------------|
| Source | Tap water |
| Conductivity* | <2000 µs/cm |
| Hardness** | <450 ppm as CaCO ₃ |
| Pressure | 0.05–0.5 MPa (7–72 psi) |
| Temperature | 5–40 °C |

Purification Water (Class III)

| Parameter | Specification |
|--------------------|---------------|
| Ionic Rejection | ≥95% |
| Bacteria Rejection | ≥99% |
| Conductivity | 1–20 µs/cm |
| Productivity Rate | 30 L/h |

High Quality Purification Water (Class II)

| Parameter | Specification |
|----------------------|---------------|
| Resistivity at 25 °C | 10 MΩ·cm |
| TOC | <30 ppb |
| Dissolved Organic | <0.1 ppm |
| Productivity Rate | 15 L/h |

Ultra-purification Water (Class I)

| Parameter | Specification |
|--------------------------|---------------|
| Resistivity at 25 °C | 18.2 MΩ·cm |
| Conductivity at 25 °C | 0.055 µs/cm |
| TOC Level*** | 1–5 ppb |
| Endotoxin (Pyrogens)**** | <0.001 EU/ml |
| Particulate (≥0.02 µm) | <1 pc/ml |
| Bacteria** | <0.1 cfu/ml |
| Rnase / Dnase** | Free |

Dispense Information

| Parameter | Specification |
|---------------------------|---------------|
| Manual Dispense Flow Rate | 1.5–2.0 L/min |
| Automatic Dispense Volume | 100–60000 ml |

Electrical Requirements

| Parameter | Specification |
|----------------------|---------------------|
| Electrical Voltage | 110 V / 220 V ± 10% |
| Electrical Frequency | 50 Hz / 60 Hz |

Packing Information

| Parameter | Specification |
|---|----------------|
| Net Weight - Main units | 35 kg |
| Net Weight - Water tank (30 L) | 7 kg |
| External Dimensions (W×D×H) - Main units | 315×525×570 mm |
| External Dimensions (W×D×H) - Water tank (30 L) | 380×380×595 mm |

- If feed water quality is poor (Conductivity > 1000 µs/cm), strengthened pretreatment module and RO-2 type is highly recommended
- When hardness of feed water is high (> 450 ppm as CaCO₃), 0.5 T water softener is recommended
- Dual wave UV module needs to be adopted. Also dependent on feed water, recommended feed TOC < 30 ppb
- Ultra-filtration module needs to be adopted. Feed water must meet the above requirements

PUREDROP - Water Purification Systems



Smart Comparison Table

| Feature | Puredrop Mini | Puredrop N | Puredrop NE/NET |
|-------------------------------|------------------------------|--------------------------------|---------------------------------|
| Water Grades | Class III, I | Class III, I | Class III, II, I |
| Ultrapure Water Quality | 18.2 MΩ·cm, TOC 1–5 ppb | 18.2 MΩ·cm, TOC 1–5 ppb | 18.2 MΩ·cm, TOC 1–5 ppb |
| High Purity (Class II) Stage | ✗ | ✗ | ✓ 10 MΩ·cm, TOC <30 ppb |
| Output Flow Rate (Manual) | 1.2–1.5 L/min | 1.5–2.0 L/min | 1.5–2.0 L/min |
| Productivity Rate (Class III) | 30 L/h | 30 L/h | 30 L/h |
| Productivity Rate (Class II) | N/A | N/A | 15 L/h |
| TOC Monitoring | ✗ | ✗ | ✓ |
| EDI Module | ✗ | ✗ | ✓ |
| Endotoxin Removal | ✗ | ✓ | ✓ |
| Particulate Removal (≥0.02μm) | ✓ | ✓ | ✓ |
| Bacteria Removal | <0.1 cfu/ml | <0.1 cfu/ml | <0.1 cfu/ml |
| Resistivity @25°C | 18.2 MΩ·cm | 18.2 MΩ·cm | 18.2 MΩ·cm |
| Conductivity @25°C | 0.055 μS/cm | 0.055 μS/cm | 0.055 μS/cm |
| Manual Dispense Volume | 1.2–1.5 L/min | 1.5–2.0 L/min | 1.5–2.0 L/min |
| Automatic Dispense Volume | 100–60000 ml | 100–60000 ml | 100–60000 ml |
| Operating Pressure | 0.05–0.5 MPa (7–72 psi) | 0.05–0.5 MPa (7–72 psi) | 0.05–0.5 MPa (7–72 psi) |
| Operating Temperature | 5–40°C | 5–40°C | 5–40°C |
| Tank Capacity | Internal 6L | External 30L | External 30L |
| Best For | HPLC Labs, Small Sample Labs | Pharma, Molecular Biology Labs | Multidisciplinary Research Labs |

PUREDROP - Water Purification Systems

Advanced Purity Technologies - EDI & TOC Monitoring EDI Module

EDI Module (Electrode ionization)

Working Principle :

- Combines ion exchange resins with ion-selective membranes and a direct electric current.
- The resin removes ions from water; the electric field moves these ions through membranes into a waste stream.
- This continuously regenerates the resin without chemical regeneration, making the process eco-friendly and consistent.

Benefits:

- Continuous ultrapure water production (up to 18.2 MΩ-cm) without downtime for regeneration.
- Eliminates chemical handling (no acids/alkalis for resin regeneration).
- Long service life and reduced operating costs.
- Stable performance for high-demand labs or industrial use.

TOC Monitoring Facility

Working Principle :






- Water passes through a UV oxidation chamber (usually 185 nm wavelength).
- The UV light oxidizes organics into CO₂, which dissolves into the water as carbonic acid.
- A conductivity sensor measures the increase in conductivity, which is directly proportional to TOC content.

Benefits:

- Real-time quality control for organics in purified water.
- Detects early contamination from bacterial by-products, leaching from components, or feed water issues.
- Critical for GLP/GMP compliance and meeting USP <643> / EP 2.2.44 standards.
- Enhances analytical accuracy for sensitive assays.



MESH-INDIA ANALYTICS PVT. LTD.

 **Head Office :** OFFICE NO.4, SR.NO. 7/17/1, RADHA PARK, NARHE, PUNE -411041
 **CIN :** U26516PN2024PTC233261  **GSTIN :** 27AASCM3860B1ZC
 **Contact :** 8149338050, 7219080093  **Email :** support@meshindia.co.in
9503578252, 9552998252